

May 22, 2023

Statement Regarding Storage Conditions and Shelf Life of the non-sterilized Cap2V8® Cap Products (includes all sizes):

The shelf life of these products is dependent upon the storage conditions. Products should be stored at normal room temperatures between 60-80° Fahrenheit in the original packaging. High-humidity locations should be avoided. Products should not be exposed to any direct sunlight and should be shielded from Ultraviolet light sources. It is best to store the product in its original sealed packaging until ready to use.

The shelf life of these SaniSure® Cap2V8® cap products is five (5) years from the date of manufacture if the storage conditions comply with those described above.